



AUTOLOGUS ADIPOSE TISSUE TRANSFER SYSTEM



# LIPOGEMS® AWARDED BEST NEW TECHNOLOGY IN SPORTS MEDICINE FOR 2016

BY ORTHOPEDICS THIS WEEK

Recognized for innovation in sports medicine for a novel invention with the potential to improve the standard of care. This exciting **NEW** technology is now available to you and your patients!



To learn more, contact your local representative or email us at [info@lipogems.us.com](mailto:info@lipogems.us.com)



**Lipogems was recognized for providing a creative and innovative solution that is cost-effective and has the potential to improve the standard of care in sports medicine.**

### **Lipogems. Your Complete Solution.**

**Process:** Closed-loop, mechanical device gently washes, rinses, and resizes autologous adipose tissue, while maintaining the vascular stromal niche

**Procedure:** Convenient point-of-care procedure for in-office applications or as an adjunct in a surgical setting

**Product:** Injectable through a small gauge needle and intended for the support and repair of soft tissue defects associated with sports injuries

The Lipogems® system is FDA cleared for indications in orthopaedic surgery.